

Ple34 Promoter PCR (pEMS1505)

MiniPromoter: Ple34
pEMS#: 1505
Expected product size (bp): 301

Reaction components	Vol/Rxn (µl)
H ₂ O	15.15
10X PCR buffer*	2.5
50 mM MgCl ₂ *	0.75
2.5 mM dNTPs**	2
10 µM Sense primer	1.25
10 µM Antisense primer	1.25
Taq Pol. (5 U/µl)*	0.1
DNA***	2
Total Volume of Rxn:	25

* Taq Polymerase set from Invitrogen (Cat no.18038-042)

** dNTPs from Invitrogen (Cat no.10297-018)


***Approximately 100 ng DNA used for samples, approximately 5 ng used for plasmid control

Samples run on a 2% agarose gel (containing SYBRsafe (Invitrogen Cat no. S33102))

Cycling conditions:	Step	Temp	Time	Note
	1	94°C	3 min	
	2	94°C	1 min	
	3	61°C	1 min	
	4	72°C	45 sec	repeat steps 2-4 34 times
	5	72°C	5 min	
	6	4°C	hold	

Primers:

Name	Sequence	T_m (°C)	Notes
oEMS2364	5' - GCGTATCACGAGGCCCTTTC - 3'	56.0	Sense primer in vector backbone
oEMS3630	5' - GGACCCCTTTGACCCATGAATACA - 3'	59.7	Anti-Sense primer for Ple34 in region "1".

Pleiades Promoter Project

1- Control plasmid DNA (pEMS1505)

2- WT mouse DNA

3- No Template

4- Knock-in ESC line

L- 100 bp ladder

expected band size = 301 bp